

For Research Use Only

SFMBT1 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60680-2-PBS



Basic Information

Catalog Number: 60680-2-PBS	GenBank Accession Number: BC014614	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 51460	CloneNo.: 1H6G4
Source: Mouse	UNIPROT ID: Q9UHI3	
Isotype: IgG1	Full Name: Scm-like with four mbt domains 1	
Immunogen Catalog Number: AG20417	Calculated MW: 866 aa, 98 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50979-1, SFMBT1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60680-1-PBS. Detection antibody: 60680-2-PBS. Standard: Ag20417. Range: 3.125-100 ng/mL.