

For Research Use Only

STIM2 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60684-1-PBS



Basic Information

Catalog Number: 60684-1-PBS	GenBank Accession Number: BC146661	Purification Method: Protein A Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 57620	CloneNo.: 3E1H10
Source: Mouse	UNIPROT ID: Q9P246	
Isotype: IgG2b	Full Name: stromal interaction molecule 2	
Immunogen Catalog Number: AG15530	Calculated MW: 754 aa, 85 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50986-1, STIM2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60684-1-PBS. Detection antibody: 60684-2-PBS. Standard: Ag15530. Range: 0.098-100 ng/mL.



Cytometric bead array standard curve of MP50986-2, STIM2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60684-1-PBS. Detection antibody: 60684-3-PBS. Standard: Ag15530. Range: 0.098-100 ng/mL.