

For Research Use Only

STIM2 Monoclonal antibody, PBS Only (Capture/Detector)

Catalog Number: 60684-2-PBS



Basic Information

Catalog Number:

60684-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG15530

GenBank Accession Number:

BC146661

GeneID (NCBI):

57620

UNIPROT ID:

Q9P246

Full Name:

stromal interaction molecule 2

Calculated MW:

754 aa, 85 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2G4E12

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50986-1, STIM2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60684-1-PBS. Detection antibody: 60684-2-PBS. Standard: Ag15530. Range: 0.098-100 ng/mL.



Cytometric bead array standard curve of MP50986-3, STIM2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60684-2-PBS. Detection antibody: 60684-3-PBS. Standard: Ag15530. Range: 0.098-100 ng/mL.