

For Research Use Only

CXCL1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60685-1-PBS



Basic Information

Catalog Number: 60685-1-PBS	GenBank Accession Number: BC011976	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 2919	CloneNo.: 1G6D8
Source: Mouse	UNIPROT ID: P09341	
Isotype: IgG2b	Full Name: chemokine (C-X-C motif) ligand 1 (melanoma growth stimulating activity, alpha)	
Immunogen Catalog Number: AG16761	Calculated MW: 8 kDa to 14 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50987-1, CXCL1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60685-1-PBS. Detection antibody: 60685-2-PBS. Standard:Ag16761. Range: 6.25-200 ng/mL