

For Research Use Only

IL20 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60691-3-PBS



Basic Information

Catalog Number:

60691-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG12447

GenBank Accession Number:

BC074949

GeneID (NCBI):

50604

UNIPROT ID:

Q9NYY1

Full Name:

interleukin 20

Calculated MW:

406 aa, 44 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

1G2G4

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50996-2, IL20 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60691-3-PBS. Detection antibody: 60691-2-PBS. Standard:Ag12447. Range: 1.563-100 ng/mL.