

For Research Use Only

Myelin basic protein Monoclonal antibody, PBS Only (Capture)

Catalog Number:60695-4-PBS



Basic Information

Catalog Number:

60695-4-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG0713

GenBank Accession Number:

BC008749

GeneID (NCBI):

4155

UNIPROT ID:

P02686

Full Name:

myelin basic protein

Calculated MW:

33 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2D8B3

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51005-3, MBP Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60695-4-PBS. Detection antibody: 60695-1-PBS. Standard:Ag0713. Range: 0.098-100 ng/mL



Cytometric bead array standard curve of MP51005-4, MBP Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60695-4-PBS. Detection antibody: 60695-3-PBS. Standard:Ag0713. Range: 0.098-100 ng/mL