

For Research Use Only

CYP4F11 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60699-1-PBS



Basic Information

Catalog Number: 60699-1-PBS	GenBank Accession Number: BC016853	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 57834	CloneNo.: 2C7D2
Source: Mouse	UNIPROT ID: Q9HBI6	
Isotype: IgG1	Full Name: cytochrome P450, family 4, subfamily F, polypeptide 11	
Immunogen Catalog Number: AG30792	Calculated MW: 524 aa, 60 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51010-1, CYP4F11 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60699-1-PBS. Detection antibody: 60699-2-PBS. Standard:Ag30792. Range: 0.195-100 ng/mL.