

For Research Use Only

CYP4F11 Monoclonal antibody, PBS Only (Detector)

Catalog Number:60699-2-PBS



Basic Information

Catalog Number:

60699-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG30792

GenBank Accession Number:

BC016853

GeneID (NCBI):

57834

UNIPROT ID:

Q9HBI6

Full Name:

cytochrome P450, family 4, subfamily

F, polypeptide 11

Calculated MW:

524 aa, 60 kDa

Purification Method:

Protein A purification

CloneNo.:

1D6B6

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51010-1, CYP4F11 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60699-1-PBS. Detection antibody: 60699-2-PBS. Standard:Ag30792. Range: 0.195-100 ng/mL.