

For Research Use Only

LTBP2 Monoclonal antibody, PBS Only (Detector)

Catalog Number:60701-4-PBS



Basic Information

Catalog Number: 60701-4-PBS	GenBank Accession Number: BC078659	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 4053	CloneNo.: 3B11G5
Source: Mouse	UNIPROT ID: Q14767	
Isotype: IgG1	Full Name: latent transforming growth factor beta binding protein 2	
Immunogen Catalog Number: AG12718	Calculated MW: 1821 aa, 195 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51013-2, LTBP2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60701-3-PBS. Detection antibody: 60701-4-PBS. Standard: Ag12718. Range: 0.391-50 ng/mL