

For Research Use Only

GBP4 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60705-4-PBS



Basic Information

Catalog Number: 60705-4-PBS	GenBank Accession Number: BC070055	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 115361	CloneNo.: 3D4C7
Source: Mouse	UNIPROT ID: Q96PP9	
Isotype: IgG1	Full Name: guanylate binding protein 4	
Immunogen Catalog Number: AG12223	Calculated MW: 640 aa, 73 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51018-3, GBP4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60705-1-PBS. Detection antibody: 60705-4-PBS. Standard:Ag12223. Range: 0.781-100 ng/mL.