

For Research Use Only

# GAP43 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60708-5-PBS



## Basic Information

|                                      |  |  |
|--------------------------------------|--|--|
| Catalog Number:<br>60708-5-PBS       | GenBank Accession Number:<br>BC007936      | Purification Method:<br>Protein G purification |
| Size:<br>1 mg/ml                     | GeneID (NCBI):<br>2596                     | CloneNo.:<br>2A2F3                             |
| Source:<br>Mouse                     | UNIPROT ID:<br>P17677                      |  |
| Isotype:<br>IgG1                     | Full Name:<br>growth associated protein 43 |  |
| Immunogen Catalog Number:<br>AG31856 | Calculated MW:<br>238 aa, 25 kDa           |  |

## Applications

Tested Applications:  
Cytometric bead array, Indirect ELISA

Species Specificity:  
human

## Background Information

### Storage

Storage:  
Store at -80°C.  
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP51022-4, GAP43 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60708-1-PBS. Detection antibody: 60708-5-PBS. Standard:Ag31856. Range: 1.563-100 ng/mL.