

For Research Use Only

RPGRIP1 Monoclonal antibody, PBS Only (Detector)

Catalog Number:60709-2-PBS



Basic Information

Catalog Number:

60709-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG31294

GenBank Accession Number:

BC039089

GeneID (NCBI):

57096

UNIPROT ID:

Q96KN7

Full Name:

retinitis pigmentosa GTPase regulator

interacting protein 1

Calculated MW:

147 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2G4H10

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51024-1, RPGRIP1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60709-1-PBS. Detection antibody: 60709-2-PBS. Standard: Ag31294. Range: 0.391-12.5 ng/mL



Cytometric bead array standard curve of MP51024-2, RPGRIP1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60709-3-PBS. Detection antibody: 60709-2-PBS. Standard: Ag31294. Range: 0.391-100 ng/mL