

For Research Use Only

# NFKBIA Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60712-1-PBS



## Basic Information

Catalog Number:

60712-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG0345

GenBank Accession Number:

BC002601

GeneID (NCBI):

4792

UNIPROT ID:

P25963

Full Name:

nuclear factor of kappa light  
polypeptide gene enhancer in B-cells  
inhibitor, alpha

Calculated MW:

36 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1E1E11

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

## Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP51027-2, NFKBIA Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66418-1-PBS. Detection antibody: 60712-1-PBS. Standard: Ag0345. Range: 0.391-100 ng/mL.



Cytometric bead array standard curve of MP51027-3, NFKBIA Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66418-2-PBS. Detection antibody: 60712-1-PBS. Standard: Ag0345. Range: 0.781-100 ng/mL.