

For Research Use Only

TNFRSF11B Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60715-2-PBS



Basic Information

Catalog Number:

60715-2-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG2139

GenBank Accession Number:

BC030155

GeneID (NCBI):

4982

UNIPROT ID:

O00300

Full Name:

tumor necrosis factor receptor
superfamily, member 11b

Calculated MW:

401 aa, 46 kDa

Purification Method:

Protein G purification

CloneNo.:

5D5B1

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51028-1, TNFRSF11B Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60715-1-PBS. Detection antibody: 60715-2-PBS. Standard:Ag2139. Range: 0.781-100 ng/mL