

For Research Use Only

SAA2 Monoclonal antibody, PBS Only (Capture)

Catalog Number:60718-1-PBS



Basic Information

Catalog Number:

60718-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG15856

GenBank Accession Number:

BC020795

GeneID (NCBI):

6289

UNIPROT ID:

P0DJ19

Full Name:

serum amyloid A2

Calculated MW:

122 aa, 14 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

2G5C2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51030-1, SAA2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60718-1-PBS. Detection antibody: 60718-2-PBS. Standard: Ag15856. Range: 0.195-100 ng/mL.