

For Research Use Only

IGHG4 Monoclonal antibody, PBS Only (Capture)

Catalog Number:60720-3-PBS



Basic Information

Catalog Number:

60720-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG9755

GenBank Accession Number:

BC025985

GeneID (NCBI):

3503

UNIPROT ID:

P01861

Full Name:

immunoglobulin heavy constant

gamma 4 (G4m marker)

Calculated MW:

473 aa, 52 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1B4B9

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51033-2, IGHG4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60720-3-PBS. Detection antibody: 60720-4-PBS. Standard: Ag9755. Range: 0.098-100 ng/mL.



Cytometric bead array standard curve of MP51033-3, IGHG4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60720-3-PBS. Detection antibody: 60720-5-PBS. Standard: Ag9755. Range: 0.098-100 ng/mL.