

For Research Use Only

Transthyretin Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60721-3-PBS



Basic Information

Catalog Number: 60721-3-PBS	GenBank Accession Number: NM_000371.4	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 7276	CloneNo.: 2D7D5
Source: Mouse	UNIPROT ID: P02766	
Isotype: IgG1	Full Name: transthyretin	
	Calculated MW: 16 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51034-2, Transthyretin Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60721-1-PBS. Detection antibody: 60721-3-PBS. Standard:Eg0931. Range: 0.098-100 ng/mL