

For Research Use Only

Tyrosinase Monoclonal antibody, PBS Only (Capture)

Catalog Number:60722-1-PBS



Basic Information

Catalog Number: 60722-1-PBS	GenBank Accession Number: BC027179	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 7299	CloneNo.: 1G1D3
Source: Mouse	UNIPROT ID: P14679	
Isotype: IgG2a	Full Name: tyrosinase (oculocutaneous albinism IA)	
Immunogen Catalog Number: AG17314	Calculated MW: 529 aa, 60 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51038-1, Tyrosinase Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60722-1-PBS. Detection antibody: 60722-2-PBS. Standard: Ag17314. Range: 0.098-12.5 ng/mL