

For Research Use Only

Tyrosinase Monoclonal antibody, PBS Only (Detector)

Catalog Number:60722-2-PBS



Basic Information

Catalog Number:

60722-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG17314

GenBank Accession Number:

BC027179

GeneID (NCBI):

7299

UNIPROT ID:

P14679

Full Name:

tyrosinase (oculocutaneous albinism IA)

Calculated MW:

529 aa, 60 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2E7H11

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51038-1, Tyrosinase Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60722-1-PBS. Detection antibody: 60722-2-PBS. Standard: Ag17314. Range: 0.098-12.5 ng/mL