

For Research Use Only

SLC20A1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60724-2-PBS



Basic Information

Catalog Number: 60724-2-PBS	GenBank Accession Number: BC019944	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 6574	CloneNo.: 2D5E3
Source: Mouse	UNIPROT ID: Q8WUM9	
Isotype: IgG1	Full Name: solute carrier family 20 (phosphate transporter), member 1	
Immunogen Catalog Number: AG28848	Calculated MW: 679 aa, 74 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51039-2, SLC20A1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60724-2-PBS. Detection antibody: 60725-2-PBS. Standard: Ag3084. Range: 0.098-100 ng/mL.