

For Research Use Only

SLC16A13 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60727-2-PBS



Basic Information

Catalog Number:

60727-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG16681

GenBank Accession Number:

BC142618

GeneID (NCBI):

201232

UNIPROT ID:

Q7RTY0

Full Name:

solute carrier family 16, member 13
(monocarboxylic acid transporter 13)

Calculated MW:

426 aa, 45 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2D8A11

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51042-1, SLC16A13 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60727-1-PBS. Detection antibody: 60727-2-PBS. Standard:Ag16681. Range: 0.781-200 ng/mL