

For Research Use Only

JUB Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60728-2-PBS



Basic Information

Catalog Number:

60728-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG29742

GenBank Accession Number:

BC007580

GeneID (NCBI):

84962

UNIPROT ID:

Q96IF1

Full Name:

jub, ajuba homolog (Xenopus laevis)

Calculated MW:

57 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

3B2G2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51044-1, JUB Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60728-1-PBS. Detection antibody: 60728-2-PBS. Standard:Ag29742. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP51044-2, JUB Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60728-3-PBS. Detection antibody: 60728-2-PBS. Standard:Ag29742. Range: 0.781-100 ng/mL.