

For Research Use Only

RABL5 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60734-5-PBS



Basic Information

Catalog Number: 60734-5-PBS	GenBank Accession Number: BC038668	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 64792	CloneNo.: 3D8A6
Source: Mouse	UNIPROT ID: Q9H7X7	
Isotype: IgG1	Full Name: RAB, member RAS oncogene family- like 5	
Immunogen Catalog Number: AG26078		

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51060-3, RABL5 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60734-5-PBS. Detection antibody: 60734-6-PBS. Standard: Ag26078. Range: 0.098-100 ng/mL.



Cytometric bead array standard curve of MP51060-4, RABL5 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60734-5-PBS. Detection antibody: 60734-7-PBS. Standard: Ag26078. Range: 0.098-100 ng/mL.