

For Research Use Only

RABL5 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60734-6-PBS



Basic Information

Catalog Number:

60734-6-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG26078

GenBank Accession Number:

BC038668

GeneID (NCBI):

64792

UNIPROT ID:

Q9H7X7

Full Name:

RAB, member RAS oncogene family-like 5

Purification Method:

Protein G purification

CloneNo.:

2F3D4

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51060-3, RABL5 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60734-5-PBS. Detection antibody: 60734-6-PBS. Standard:Ag26078. Range: 0.098-100 ng/mL.