

For Research Use Only

FUT4 Monoclonal antibody, PBS Only (Capture)

Catalog Number:60736-1-PBS



Basic Information

Catalog Number:

60736-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG17695

GenBank Accession Number:

BC136374

GeneID (NCBI):

2526

UNIPROT ID:

P22083

Full Name:

fucosyltransferase 4 (alpha (1,3)
fucosyltransferase, myeloid-specific)

Calculated MW:

530 aa, 59 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1H3E4

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51063-1, FUT4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60736-1-PBS. Detection antibody: 60736-2-PBS. Standard:Ag17695. Range: 0.098-100 ng/mL



Cytometric bead array standard curve of MP51063-2, FUT4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60736-1-PBS. Detection antibody: 60736-3-PBS. Standard:Ag17695. Range: 0.098-100 ng/mL