

For Research Use Only

# LRIG1 Monoclonal antibody, PBS Only (Detector)

Catalog Number:60738-4-PBS



## Basic Information

Catalog Number:

60738-4-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG17880

GenBank Accession Number:

BC071561

GeneID (NCBI):

26018

UNIPROT ID:

Q96JA1

Full Name:

leucine-rich repeats and immunoglobulin-like domains 1

Calculated MW:

1093 aa, 119 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

1D6C4

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP51067-2, LRIg1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60738-3-PBS. Detection antibody: 60738-4-PBS. Standard:Ag17880. Range: 0.781-100 ng/mL.