

For Research Use Only

LRIG1 Monoclonal antibody, PBS Only (Detector)

Catalog Number:60738-5-PBS



Basic Information

Catalog Number: 60738-5-PBS	GenBank Accession Number: BC071561	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 26018	CloneNo.: 2H8E7
Source: Mouse	UNIPROT ID: Q96JA1	
Isotype: IgG1	Full Name: leucine-rich repeats and immunoglobulin-like domains 1	
Immunogen Catalog Number: AG17880	Calculated MW: 1093 aa, 119 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51067-3, LRIg1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60738-3-PBS. Detection antibody: 60738-5-PBS. Standard:Ag17880. Range: 0.781-100 ng/mL.