

For Research Use Only

WDR78 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60742-4-PBS



Basic Information

Catalog Number: 60742-4-PBS	GenBank Accession Number: BC126379	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 79819	CloneNo.: 2C4A10
Source: Mouse	UNIPROT ID: Q5VTH9	
Isotype: IgG2a	Full Name: WD repeat domain 78	
Immunogen Catalog Number: AG27018		

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51072-3, WDR78 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60742-4-PBS. Detection antibody: 60742-5-PBS. Standard:Ag27018. Range: 0.195-100 ng/mL.