

For Research Use Only

# IFT52 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60743-3-PBS



## Basic Information

Catalog Number:

60743-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG11613

GenBank Accession Number:

BC039831

GeneID (NCBI):

51098

UNIPROT ID:

Q9Y366

Full Name:

intraflagellar transport 52 homolog (Chlamydomonas)

Calculated MW:

437 aa, 50 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

4G5G5

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP51076-2, IFT52 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60743-1-PBS. Detection antibody: 60743-3-PBS. Standard:Ag11613. Range: 0.391-100 ng/mL.