

For Research Use Only

# SGK1 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60747-1-PBS



## Basic Information

|  |   |   |
|--|---|---|
| <b>Catalog Number:</b><br>60747-1-PBS      | <b>GenBank Accession Number:</b><br>BC001263                    | <b>Purification Method:</b><br>Protein A purification |
| <b>Size:</b><br>1 mg/ml                    | <b>GeneID (NCBI):</b><br>6446                                   | <b>CloneNo.:</b><br>1H7G9                             |
| <b>Source:</b><br>Mouse                    | <b>UNIPROT ID:</b><br>O00141                                    |   |
| <b>Isotype:</b><br>IgG3                    | <b>Full Name:</b><br>serum/glucocorticoid regulated<br>kinase 1 |   |
| <b>Immunogen Catalog Number:</b><br>AG2411 | <b>Calculated MW:</b><br>431 aa, 49 kDa                         |   |

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP51079-1, SGK1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60746-1-PBS. Detection antibody: 60747-1-PBS. Standard:Ag2411. Range: 0.391-12.5 ng/mL



Cytometric bead array standard curve of MP51079-2, SGK1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60747-2-PBS. Detection antibody: 60747-1-PBS. Standard:Ag2411. Range: 0.391-100 ng/mL