

For Research Use Only

CCL22/MDC Monoclonal antibody, PBS Only (Capture/Detector)

Catalog Number:60754-3-PBS



Basic Information

Catalog Number:

60754-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG4275

GenBank Accession Number:

BC027952

GeneID (NCBI):

6367

UNIPROT ID:

O00626

Full Name:

chemokine (C-C motif) ligand 22

Calculated MW:

93 aa, 11 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

2B8G7

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51202-2, CCL22/MDC Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60754-3-PBS. Detection antibody: 60754-2-PBS. Standard:Ag4275. Range: 0.781-50 ng/mL.



Cytometric bead array standard curve of MP51202-4, CCL22/MDC Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60754-5-PBS. Detection antibody: 60754-3-PBS. Standard:Ag4275. Range: 0.781-100 ng/mL.