

For Research Use Only

CCL22/MDC Monoclonal antibody, PBS Only (Capture)

Catalog Number:60754-4-PBS



Basic Information

Catalog Number: 60754-4-PBS	GenBank Accession Number: BC027952	Purification Method: Protein A Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 6367	CloneNo.: 1E12E2
Source: Mouse	UNIPROT ID: O00626	
Isotype: IgG2b	Full Name: chemokine (C-C motif) ligand 22	
Immunogen Catalog Number: AG4275	Calculated MW: 93 aa, 11 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51202-3, CCL22/MDC Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60754-4-PBS. Detection antibody: 60754-2-PBS. Standard:Ag4275. Range: 0.781-100 ng/mL