

For Research Use Only

Cytokeratin 17 Monoclonal antibody, PBS Only (Detector)

Catalog Number:60771-6-PBS



Basic Information

Catalog Number: 60771-6-PBS	GenBank Accession Number: BC000159	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 3872	CloneNo.: 1F9E3
Source: Mouse	UNIPROT ID: Q04695	
Isotype: IgG1	Full Name: keratin 17	
Immunogen Catalog Number: AG17657	Calculated MW: 432 aa, 48 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51099-3, Cytokeratin 17 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60771-5-PBS. Detection antibody: 60771-6-PBS. Standard:Ag17657. Range: 1.563-200 ng/mL.



Cytometric bead array standard curve of MP51099-4, Cytokeratin 17 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60771-3-PBS. Detection antibody: 60771-6-PBS. Standard:Ag17657. Range: 1.563-200 ng/mL.