

For Research Use Only

LAMC3-Specific Monoclonal antibody, PBS Only (Capture)

Catalog Number:60773-1-PBS



Basic Information

Catalog Number:

60773-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG19466

GenBank Accession Number:

NM_006059

GeneID (NCBI):

10319

UNIPROT ID:

Q9Y6N6

Full Name:

laminin, gamma 3

Calculated MW:

1575 aa, 171 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

2D7B12

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51059-3, LAMC3-Specific Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60773-1-PBS. Detection antibody: 60773-2-PBS. Standard: Ag19466. Range: 0.781-12.5 ng/mL