

For Research Use Only

P21;CDKN1A Monoclonal antibody, PBS Only (Capture)

Catalog Number:60775-3-PBS



Basic Information

Catalog Number:

60775-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG18100

GenBank Accession Number:

BC000275

GeneID (NCBI):

1026

UNIPROT ID:

P38936

Full Name:

cyclin-dependent kinase inhibitor 1A
(p21, Cip1)

Calculated MW:

18 kDa

Purification Method:

Protein G purification

CloneNo.:

3C1C10

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51104-3, P21 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60775-3-PBS. Detection antibody: 60775-4-PBS. Standard: Ag0308. Range: 1.563-25 ng/mL.