

For Research Use Only

CTNND2 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60780-1-PBS



Basic Information

Catalog Number:

60780-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG28041

GenBank Accession Number:

NM_001288715

GeneID (NCBI):

1501

UNIPROT ID:

Q9UQB3

Full Name:

catenin (cadherin-associated protein),
delta 2 (neural plakophilin-related
arm-repeat protein)

Calculated MW:

133 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1E5G1

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51111-1, CTNND2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60780-1-PBS. Detection antibody: 60780-2-PBS. Standard:Ag28041. Range: 0.781-50 ng/mL