

For Research Use Only

# Claudin 18 Monoclonal antibody, PBS Only (Detector)

Catalog Number:60782-2-PBS



## Basic Information

Catalog Number: 60782-2-PBS	GenBank Accession Number: NM_001002026	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 51208	CloneNo.: 1A12E2
Source: Mouse	UNIPROT ID: P56856-2	
Isotype: IgG1	Full Name: claudin 18	
Immunogen Catalog Number: AG34049	Calculated MW: 28 kDa	

## Applications

Tested Applications:  
Cytometric bead array, Indirect ELISA

Species Specificity:  
human

## Background Information

### Storage

Storage:  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP51118-1, Claudin 18 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60782-1-PBS. Detection antibody: 60782-2-PBS. Standard: Ag34049. Range: 6.25-100 ng/mL.



Cytometric bead array standard curve of MP51118-2, Claudin 18 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60782-3-PBS. Detection antibody: 60782-2-PBS. Standard: Ag34049. Range: 3.125-100 ng/mL.