

For Research Use Only

# CSF2RB Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60786-2-PBS



## Basic Information

Catalog Number:

60786-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG28159

GenBank Accession Number:

BC070085

GeneID (NCBI):

1439

UNIPROT ID:

P32927

Full Name:

colony stimulating factor 2 receptor,  
beta, low-affinity (granulocyte-  
macrophage)

Calculated MW:

889 aa, 96 kDa

Purification Method:

Protein A purification

CloneNo.:

1G2A2

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP51124-1, CSF2RB Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60786-1-PBS. Detection antibody: 60786-2-PBS. Standard:Ag28159. Range: 0.781-100 ng/mL.