

For Research Use Only

MARVELD3 Monoclonal antibody, PBS Only (Detector)

Catalog Number:60788-2-PBS



Basic Information

Catalog Number: 60788-2-PBS	GenBank Accession Number: BC013376	Purification Method: Protein A Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 91862	CloneNo.: 1F6E7
Source: Mouse	UNIPROT ID: Q96A59	
Isotype: IgG2b	Full Name: MARVEL domain containing 3	
Immunogen Catalog Number: AG33294	Calculated MW: 401 aa, 45 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51126-1, MARVELD3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60788-1-PBS. Detection antibody: 60788-2-PBS. Standard: Ag33294. Range: 3.125-100 ng/mL.