

For Research Use Only

# LAMA1 Monoclonal antibody, PBS Only (Capture)

Catalog Number:60790-3-PBS



## Basic Information

Catalog Number: 60790-3-PBS	GenBank Accession Number: NM_005559	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 284217	CloneNo.: 1E9A3
Source: Mouse	UNIPROT ID: P25391	
Isotype: IgG2a	Full Name: laminin, alpha 1	
Immunogen Catalog Number: AG29326	Calculated MW: 337 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP51131-2, LAMA1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60790-3-PBS. Detection antibody: 60790-4-PBS. Standard:Ag29326. Range: 0.098-100 ng/mL.