

For Research Use Only

Plakophilin 2 Monoclonal antibody, PBS Only (Capture)

Catalog Number:60795-1-PBS



Basic Information

Catalog Number: 60795-1-PBS	GenBank Accession Number: BC094762	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 5318	CloneNo.: 3C8E11
Source: Mouse	UNIPROT ID: Q99959	
Isotype: IgG1	Full Name: plakophilin 2	
Immunogen Catalog Number: AG30401	Calculated MW: 837 aa, 93 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51142-1, PKP2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60795-1-PBS. Detection antibody: 60795-2-PBS. Standard:Ag30401. Range: 0.195-200 ng/mL.