

For Research Use Only

# TTBK2 Monoclonal antibody, PBS Only (Detector)

Catalog Number:60796-4-PBS



## Basic Information

Catalog Number: 60796-4-PBS	GenBank Accession Number: BC071556	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 146057	CloneNo.: 3C3A11
Source: Mouse	UNIPROT ID: Q6IQ55	
Isotype: IgG1	Full Name: tau tubulin kinase 2	
Immunogen Catalog Number: AG30470	Calculated MW: 137 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP51144-3, TTBK2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60796-1-PBS. Detection antibody: 60796-4-PBS. Standard: Ag30470. Range: 0.781-100 ng/mL.