

For Research Use Only

CD247 Monoclonal antibody, PBS Only (Capture/Detector)

Catalog Number: 60797-4-PBS



Basic Information

Catalog Number: 60797-4-PBS	GenBank Accession Number: BC025703	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 919	CloneNo.: 3C8F2
Source: Mouse	UNIPROT ID: P20963	
Isotype: IgG1	Full Name: CD247 molecule	
Immunogen Catalog Number: AG33870	Calculated MW: 16 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51145-2, CD247 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60797-3-PBS. Detection antibody: 60797-4-PBS. Standard:Ag33870. Range: 0.391-100 ng/mL.



Cytometric bead array standard curve of MP51145-3, CD247 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60797-4-PBS. Detection antibody: 60797-2-PBS. Standard:Ag33870. Range: 0.391-100 ng/mL.