

For Research Use Only

KRT10 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60798-1-PBS



Basic Information

Catalog Number: 60798-1-PBS	GenBank Accession Number: NM_000421	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 3858	CloneNo.: 3G6E8
Source: Mouse	UNIPROT ID: P13645	
Isotype: IgG1	Full Name: keratin 10	
Immunogen Catalog Number: AG32350	Calculated MW: 59kd	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51147-1, Cytokeratin 10 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60798-1-PBS. Detection antibody: 60798-2-PBS. Standard:Ag32350. Range: 0.195-100 ng/mL



Cytometric bead array standard curve of MP51147-2, Cytokeratin 10 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60798-1-PBS. Detection antibody: 60798-3-PBS. Standard:Ag32350. Range: 0.391-100 ng/mL