

For Research Use Only

B3GALNT2 Monoclonal antibody, PBS Only (Detector)

Catalog Number:60799-2-PBS



Basic Information

Catalog Number: 60799-2-PBS	GenBank Accession Number: BC029564	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 148789	CloneNo.: 2F2E3
Source: Mouse	UNIPROT ID: Q8NCRO	
Isotype: IgG1	Full Name: beta-1,3-N-acetylgalactosaminyltransferase 2	
Immunogen Catalog Number: AG31400	Calculated MW: 500 aa, 57 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51148-1, B3GALNT2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60799-1-PBS. Detection antibody: 60799-2-PBS. Standard: Ag31400. Range: 12.5-200 ng/mL