

For Research Use Only

EPB41 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60801-2-PBS



Basic Information

Catalog Number: 60801-2-PBS	GenBank Accession Number: BC039079	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 2035	CloneNo.: 1A12C9
Source: Mouse	UNIPROT ID: P11171	
Isotype: IgG2a	Full Name: erythrocyte membrane protein band 4.1 (elliptocytosis 1, RH-linked)	
Immunogen Catalog Number: AG4365	Calculated MW: 720 aa, 81 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51150-1, EPB41 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60801-1-PBS. Detection antibody: 60801-2-PBS. Standard: Ag4365. Range: 3.125-100 ng/mL.