

For Research Use Only

LEF1 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60803-2-PBS



Basic Information

Catalog Number:

60803-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG29629

GenBank Accession Number:

BC050632

GeneID (NCBI):

51176

ENSEMBL Gene ID:

ENSG00000138795

UNIPROT ID:

Q9UJU2

Full Name:

Lymphoid enhancer-binding factor 1

Calculated MW:

37 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1F3H5

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51153-1, LEF1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60803-1-PBS. Detection antibody: 60803-2-PBS. Standard: Ag29629. Range: 1.563-25 ng/mL