

For Research Use Only

AOC2 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60806-1-PBS



Basic Information

Catalog Number: 60806-1-PBS	GenBank Accession Number: BC142641	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 314	CloneNo.: 3F3G8
Source: Mouse	UNIPROT ID: O75106	
Isotype: IgG2b	Full Name: amine oxidase, copper containing 2 (retina-specific)	
Immunogen Catalog Number: AG26965	Calculated MW: 756 aa, 84 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51161-1, AOC2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60806-1-PBS. Detection antibody: 60806-2-PBS. Standard: Ag26965. Range: 0.391-50 ng/mL