

For Research Use Only

# JUP Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60807-1-PBS



## Basic Information

Catalog Number:

60807-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgM

Immunogen Catalog Number:

AG17598

GenBank Accession Number:

BC011865

GeneID (NCBI):

3728

UNIPROT ID:

P14923

Full Name:

junction plakoglobin

Calculated MW:

745 aa, 82 kDa

Purification Method:

Thiophilic affinity chromatograph

CloneNo.:

1B8E11

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP51165-1, JUP Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60807-1-PBS. Detection antibody: 66445-2-PBS. Standard:Ag1627. Range: 1.563-100 ng/mL.



Cytometric bead array standard curve of MP51165-2, JUP Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60807-1-PBS. Detection antibody: 66445-3-PBS. Standard:Ag1627. Range: 3.125-100 ng/mL.