

For Research Use Only

APLP1 Monoclonal antibody, PBS Only (Detector)

Catalog Number:60819-2-PBS



Basic Information

| | | |
|---|---|---|
| Catalog Number: 60819-2-PBS | GenBank Accession Number: BC012889 | Purification Method: Protein G purification |
| Size: 1 mg/ml | GeneID (NCBI): 333 | CloneNo.: 1C3F10 |
| Source: Mouse | UNIPROT ID: P51693 | |
| Isotype: IgG1 | Full Name: amyloid beta (A4) precursor-like protein 1 | |
| Immunogen Catalog Number: AG28632 | Calculated MW: 650 aa, 72 kDa | |

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51193-1, APLP1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60819-1-PBS. Detection antibody: 60819-2-PBS. Standard: Ag28632. Range: 1.566-100 ng/mL.